



RETRACTOR MAKER

PRODUCT BROCHURE OF

Judd-Masson Self Retaining Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Judd-Masson Self Retaining Retractor

SKU RM-JMSR-399-1

\$399 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Heavy-Duty Self-Retaining Base Frame	Multi-point chassis engineered for wide 3D visibility
2	Fixed Flat Lateral Blades	69 mm deep outward-facing panels with 90-degree protective fenestrations
1	Adjustable Central Blade Component	Solid flat blade configured to safely contort to internal organ structures
1	Dual Manual Ratchet Lock Assembly	Dual mechanical locking clips managing expansion settings for both blade systems
2	Ergonomic Retracting Forceps Beams	Lightweight structured handles built with smooth tissue-protective edges
1	Protective Storage / Box	Standard transit packaging

PRODUCT DESCRIPTION

The Judd-Masson Self Retaining Retractor is an internationally acclaimed urological and pelvic surgical retractor system engineered specifically to provide safe, multi-directional soft-tissue and visceral displacement during open suprapubic urological procedures. Designed with an ultra-rigid multi-point chassis, this instrument holds the abdominal wall, rectus muscles, and bladder margins open under locked mechanical tension, eliminating assistant hand fatigue while preserving a wide, static visual corridor into the bladder lumen.

Built to clinical framework specifications, this Judd-Masson Self Retaining Retractor complete set (\$399) features a heavy-duty self-retaining base frame engineered for wide 3D visibility inside pelvic and abdominal spaces. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the system features a 68 mm

standard operative structural depth paired with a lateral blade depth profile of 69 mm deep outward-facing lateral tissue panels with 90-degree protective fenestrations. Its adjustable solid flat central blade safely contorts to deep internal organ structures, while the dual manual ratchet lock assembly manages expansion settings for both blade systems up to a 120 mm maximum horizontal track limit.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Judd-Masson Self Retaining Retractor to withstand repetitive hospital autoclave cycles, providing urology operating theaters, pelvic reconstructive centers, and general surgery suites worldwide with exceptional mechanical stability, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)