



RETRACTOR MAKER

PRODUCT BROCHURE OF

Doyen Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Doyen Retractor

SKU RM-DR-22-1

\$99 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
30 x 30 mm blade — 23 cm (9½")	RM-DR-22-1-1	\$89
32 x 40 mm blade — 23 cm (9½")	RM-DR-22-1-2	\$94
40 x 50 mm blade — 23 cm (9½")	RM-DR-22-1-3	\$99
40 x 60 mm blade — 23 cm (9½")	RM-DR-22-1-4	\$104
45 x 75 mm blade — 23 cm (9½")	RM-DR-22-1-5	\$109
65 x 85 mm blade — 23 cm (9½")	RM-DR-22-1-6	\$114

PRODUCT DESCRIPTION

The Doyen Retractor is an internationally recognized general surgical instrument engineered specifically to provide controlled, atraumatic tissue displacement during open abdominal and pelvic operations. Designed with a wide, curved blade profile, this instrument distributes manual counter-traction across a broad anatomical surface, preventing localized pressure necrosis on retracted organ walls and muscle planes.

Built for optimal operational control, this Doyen Retractor model features an overall structural reach of 22 cm (\$99). Its specialized wide blade profile is contoured to optimize anatomical exposure while minimizing localized tissue pressure, allowing surgical assistants to displace the abdominal wall, visceral masses, and bladder smoothly. Its ergonomic handle incorporates a contoured finger ring that delivers non-slip manual control and superior pulling leverage throughout extended surgical cases.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Doyen Retractor to strict clinical standards, delivering high tensile

strength, exceptional corrosion resistance, and complete autoclave durability for hospital operating theaters worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)