



RETRACTOR MAKER

PRODUCT BROCHURE OF

Pelvic Retractor Blunt

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Pelvic Retractor Blunt

SKU RM-PRB-105-1

\$89 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Pelvic Retractor Blunt	10 ½" (26.67 cm) length, 1" wide solid blunt blade	\$89
Protective Sleeve / Box	Individual transit packaging	\$89

PRODUCT DESCRIPTION

The Pelvic Retractor Blunt is a premier deep-cavity surgical retractor engineered specifically to provide safe, non-traumatic displacement of pelvic viscera, bladder, and deep muscular walls during major pelvic and abdominal operations. Designed with a broad, curved blade terminating in a soft blunt tip, this instrument isolates deep pelvic organs cleanly away from the operative field without causing localized contusions or tissue lacerations.

Built to withstand frequent sterilization protocols while remaining highly resistant to corrosion, this Pelvic Retractor Blunt configuration features an overall tool length of 10 ½ inches (26.67 cm, \$89) that provides excellent reach into the true pelvis. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the working end incorporates a 1-inch (2.54 cm) wide solid blunt blade layout contoured to cradle pelvic tissues safely under continuous manual counter-traction. The ergonomic handle ensures operating assistants maintain fatigue-free manual leverage throughout extended pelvic reconstructive procedures.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Pelvic Retractor Blunt to withstand heavy mechanical loads and repeated hospital autoclave cycles, providing gynecological, urological, and orthopedic pelvic trauma suites worldwide with unyielding structural rigidity, high corrosion resistance, and complete clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)