



RETRACTOR MAKER

PRODUCT BROCHURE OF

Cooley Atrial Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Cooley Atrial Retractor

SKU RM-CAR-1

\$64 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Small Model	RM-CAR-1-1	\$64
Medium Model	RM-CAR-1-2	\$69
Large Model	RM-CAR-1-3	\$74

PRODUCT DESCRIPTION

The Cooley Atrial Retractor is a premier cardiac surgical instrument engineered specifically to provide clear, atraumatic visualization of the left atrium and mitral valve annulus during open heart procedures. Designed with an anatomically curved distal paddle, this instrument allows surgeons to elevate the atrial wall and interatrial septum cleanly, creating an unhindered line of sight for valve inspection and repair.

Built for optimal deep-cavity reach, this Cooley Atrial Retractor model features a uniform 24 cm overall structural length that keeps the assistant's hands safely away from the primary operative field. The instrument is manufactured in three selectable blade sizes: Small (\$64), Medium (\$69), and Large (\$74), allowing surgical teams to match patient-specific atrial chamber dimensions. Its slender shaft and balanced ergonomic flat handle ensure fatigue-free manual counter-traction during delicate intra-cardiac suturing.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Cooley Atrial Retractor to withstand repetitive hospital autoclave cycles, providing cardiothoracic surgery centers worldwide with superior corrosion resistance, unyielding structural rigidity, and dependable durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)