



RETRACTOR MAKER

PRODUCT BROCHURE OF

Cloward Cervical/Lamina Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Cloward Cervical/Lamina Retractor

SKU RM-CCLR-1

\$99 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Cross-Serrated Jaws	RM-CCLR-1-1	\$99
With Teeth Jaws	RM-CCLR-1-2	\$99

PRODUCT DESCRIPTION

The Cloward Cervical/Lamina Retractor is a foundational spinal surgical instrument engineered specifically to provide secure, hands-free exposure of the cervical lamina, facet joints, and paraspinal musculature during posterior spine procedures. Designed with an optimal balance of reach and mechanical leverage, this instrument eliminates assistant fatigue while maintaining an unyielding operative corridor.

Built for precise spinal navigation, this Cloward Cervical/Lamina Retractor configuration features an overall length of 16 cm (\$99) that provides adequate depth while preserving high manual control in cervical corridors. Surgeons can configure the instrument with either Cross-Serrated Jaws engineered to deliver a secure, distributed grip on soft tissue structures, or Toothed Jaws optimized to provide maximum mechanical traction and non-slip purchase on dense muscular or fascial layers. The integrated finger-ring ratchet catch maintains constant retraction tension under heavy muscle resistance.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Cloward Cervical/Lamina Retractor to withstand repetitive hospital sterilization cycles, providing spine and neurosurgical operating theaters worldwide with exceptional durability and clinical stability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)