



RETRACTOR MAKER

PRODUCT BROCHURE OF

Wylie Renal Vein Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Wylie Renal Vein Retractor

SKU RM-WRVR-1

\$89 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
11 ¾" (29.85 cm) — Smooth curved vascular blade	RM-WRVR-1-1	\$89
12.0" (30.48 cm) — Smooth curved vascular blade	RM-WRVR-1-2	\$94
12 ½" (31.75 cm) — Smooth curved vascular blade	RM-WRVR-1-3	\$99
12 ¾" (32.39 cm) — Smooth curved vascular blade	RM-WRVR-1-4	\$104
13.0" (33.02 cm) — Smooth curved vascular blade	RM-WRVR-1-5	\$109

PRODUCT DESCRIPTION

The Wylie Renal Vein Retractor is a specialized vascular and urological surgical retractor engineered specifically to provide safe, non-traumatic displacement and mobilization of the renal vein, inferior vena cava, and adjacent retroperitoneal vessels. Designed with an elongated curved blade terminating in a soft retaining contour, this instrument allows vascular surgeons to elevate large venous channels gently away from the abdominal aorta without causing endothelial shearing, tearing, or avulsion.

Built to withstand frequent autoclave sterilization cycles while maintaining high structural resilience, this Wylie Renal Vein Retractor series features an adaptive sizing matrix across five selectable overall lengths: 11 ¾ inches (29.85 cm, \$89), 12 inches (30.48 cm, \$94), 12 ½ inches (31.75 cm, \$99), 12 ¾ inches (32.39 cm, \$104), and 13 inches (33.02 cm, \$109). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic),

the instrument gives surgical assistants superior reaching depth into the retroperitoneum while keeping manual handling clear of primary vascular clamps. Its ergonomic flat handle provides fatigue-free manual leverage.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Wylie Renal Vein Retractor to uphold strict clinical standards under intense surgical lighting, providing cardiovascular suites, urological operating theaters, and renal transplant centers worldwide with exceptional corrosion resistance, structural rigidity, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)