



RETRACTOR MAKER

PRODUCT BROCHURE OF

Kelly Retractor-Lamb Handle

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Kelly Retractor-Lamb Handle

SKU RM-KLR-LH-1

\$94 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	51 mm long × 63 mm wide	Overall Length: 24 cm · Primary Clinical Role: Standard abdominal wall and intermediate muscular retraction · Price (USD): \$94.00
1	63 mm long × 76 mm wide	Overall Length: 24 cm · Primary Clinical Role: Deep pelvic floor and laparotomy visceral displacement · Price (USD): \$94.00
1	76 mm long × 89 mm wide	Overall Length: 24 cm · Primary Clinical Role: Maximum-width abdominal wall and bariatric visceral exposure · Price (USD): \$94.00
1	Protective Sleeve / Box	Overall Length: Standard Transit Packaging · Primary Clinical Role: — · Price (USD): Included

PRODUCT DESCRIPTION

The Kelly Retractor-Lamb Handle is a specialized general and gynecological surgical instrument engineered specifically to provide powerful, controlled tissue displacement during major open abdominal and pelvic operations. Designed with an ergonomic Lamb handle profile, this instrument allows surgical assistants to apply heavy manual counter-traction with optimal palm support, preventing grip slippage and hand fatigue during deep wound exposure.

Built to ensure controlled tissue displacement across variable patient anatomies, this Kelly Retractor-Lamb Handle features a uniform overall instrument length of 24 cm (\$94). The retractor is offered across three selectable blade footprints: 51 mm long × 63 mm wide, 63 mm long × 76 mm wide, and an expansive 76 mm long × 89 mm wide profile. The right-angled blade incorporates a smooth, rounded retaining lip that holds

viscera and dense muscle planes securely without puncturing organ serosa.

Crafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Kelly Retractor-Lamb Handle to withstand heavy mechanical loads and repeated hospital autoclave cycles, providing operating rooms worldwide with exceptional durability, corrosion resistance, and complete clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)