



RETRACTOR MAKER

PRODUCT BROCHURE OF

Coronal Brow Lift Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Coronal Brow Lift Retractor

SKU RM-CBLR-1

\$49 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Coronal Brow Lift Retractor	7 ¼" (18.4 cm) length, 70.0 mm rake head, 7 sharp prongs	\$49
Protective Sleeve / Box	Individual transit packaging	\$49

PRODUCT DESCRIPTION

The Coronal Brow Lift Retractor is a specialized plastic and facial reconstructive surgical instrument engineered specifically to deliver powerful, non-slip tissue traction during open coronal forehead and brow lift procedures. Designed with an expansive 70.0 mm wide rake head, this instrument distributes manual counter-traction across a broad galeal surface area, preventing focal tearing while maximizing operative visualization of the supraorbital rim.

Built for heavy leverage without slippage, this Coronal Brow Lift Retractor configuration features an overall length of 7 ¼ inches (18.4 cm, \$49). Its working end is equipped with 7 sharp offset prongs specifically engineered to grip the dense fibrous scalp, galea, and frontalis musculature firmly during continuous manual pulling. The angled shank keeps the assistant's hand outside the line of sight, while the ergonomic flat handle provides fatigue-free manual leverage throughout deep subperiosteal tunneling.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Coronal Brow Lift Retractor to withstand rigorous hospital autoclave cycles, delivering high tensile strength, exceptional corrosion resistance, and complete clinical durability for cosmetic and facial plastic surgery centers worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)