



RETRACTOR MAKER

PRODUCT BROCHURE OF

Converse Nasal Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Converse Nasal Retractor

SKU RM-CNR-45-1

\$69 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Without Fiber Optic Light	RM-CNR-45-1-1	\$69
With Fiber Optic Light	RM-CNR-45-1-2	\$209

PRODUCT DESCRIPTION

The Converse Nasal Retractor is an industry-standard rhinoplasty surgical retractor engineered specifically to provide steady, atraumatic elevation of the nasal skin envelope, cartilaginous vault, and dorsal periosteum. Designed with an optimal downward angle, this instrument allows surgeons to maintain clear operative access to the bony and cartilaginous dorsum while protecting the overlying skin flap from focal bruising.

Built for procedural adaptability, this Converse Nasal Retractor model features an overall length of 4 ½ inches (11.4 cm) paired with an 11.0 mm wide by 44.0 mm long tapered blade profile. The instrument is available in two functional options: a standard non-illuminated model (\$69) and a premium illuminated configuration (\$209) incorporating an integrated fiber-optic light bundle that projects bright, shadowless illumination directly into the deep surgical pocket. The light bundle includes a universal fitting compatible with standard OR light sources (Storz, Olympus, Wolf, ACMI).

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Converse Nasal Retractor to withstand repetitive hospital sterilization cycles, providing plastic and ENT surgical centers with exceptional optical clarity, corrosion resistance, and clinical longevity.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)