



RETRACTOR MAKER

PRODUCT BROCHURE OF

Adjustable Mayfield Adaptor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Adjustable Mayfield Adaptor

SKU RM-AMA-1

\$1,199 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Adjustable Mayfield Adaptor	Provides a secure mounting platform for Mayfield-style cranial fixation systems on operating tables without built-in receptacle holes
2	Adjustable End Brackets	Secure the adaptor to operating table side rails and accommodate different table widths
Multiple	Mounting Holes	Allow lateral adjustment and flexible positioning of the skull clamp base unit

PRODUCT DESCRIPTION

The Adjustable Mayfield Adaptor is an essential surgical table mounting bridge engineered to enable rigid cranial fixation on operating tables that do not feature built-in skull clamp receptacle sockets. By establishing an ultra-secure mechanical connection across the table's accessory side rails, this adaptor ensures that cranial fixation frames maintain absolute rigidity without play or slippage.

Equipped with quick-adjusting telescoping end brackets, the Adjustable Mayfield Adaptor spans table widths from 21.375 inches to 24.5 inches (54.3 cm to 62.2 cm) with positive-locking side rail clamps. Its crossbar platform incorporates multiple precision-bored mounting holes, allowing the surgical team to adjust the skull clamp base unit laterally to accommodate off-center patient positioning, park-bench setups, or lateral decubitus alignments.

CNC-machined from high-tensile, aircraft-grade metal and aluminum alloys with hard-coat anodized protection, Retractor Maker supplies the Adjustable Mayfield Adaptor non-sterile and fully autoclavable inside an export-grade protective box with high-density custom foam inserts.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)