



RETRACTOR MAKER

PRODUCT BROCHURE OF

Neivert Knife and Guide Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Neivert Knife and Guide Retractor

SKU RM-NKGR-1

\$39 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Small — 6 3/4 in	RM-NKGR-1-1	\$39
Standard — 7 3/4 in	RM-NKGR-1-2	\$44

PRODUCT DESCRIPTION

The Neivert Knife and Guide Retractor is an indispensable rhinoplasty and facial plastic surgical instrument engineered specifically to provide precise soft-tissue retraction and controlled osteotomy knife guidance during endonasal interventions. Designed with a dual-ended layout, this instrument allows surgeons to elevate delicate alar cartilage flaps with one end and seamlessly guide intercartilaginous knife passes using the opposing channel without switching tools.

Built for complete procedural adaptability, this Neivert Knife and Guide Retractor holds a fixed dual-ended tip layout available across two selectable overall tool lengths: 6 ¾ inches (17.15 cm, \$39) and 7 ¾ inches (19.69 cm, \$44). Working Tip 1 features a blunted 15 mm skin retractor equipped with a rounded lip to prevent accidental cutaneous tearing. Working Tip 2 features a 4 mm closed knife guide channel layout engineered to stabilize scalpel blades during precise cartilage division. The central flat handle delivers optimal ergonomic balance and tactile feedback.

Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Neivert Knife and Guide Retractor to withstand repetitive hospital autoclave cycles, providing plastic, reconstructive, and ENT surgical suites worldwide with exceptional corrosion resistance, structural rigidity, and clinical longevity.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)