



RETRACTOR MAKER

PRODUCT BROCHURE OF

Love Nerve Root Retractor, Bayonet

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Love Nerve Root Retractor, Bayonet

SKU RM-LNRR-BAY-21-1

\$29 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Love Nerve Root Retractor (Bayonet)	Dimensions / Details: 21 cm overall length, bayonet shaft, cupped nerve blade · Quantity Included: 1	\$29
Protective Sleeve / Box	Dimensions / Details: Individual transit packaging · Quantity Included: 1	\$29

PRODUCT DESCRIPTION

The Love Nerve Root Retractor is a specialized microsurgical spine instrument engineered specifically to provide safe, non-traumatic nerve root displacement while keeping the surgeon's hand and handle completely outside the microscopic optical axis. Designed with an offset bayonet shaft geometry, this instrument drops the manual handling corridor below the visual line of sight, allowing continuous, unobstructed visualization of the intervertebral disc space.

Optimized for complex spine environments, this Love Nerve Root Retractor configuration features an overall tool length of 21 cm (\$29). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument incorporates a bayonet configuration profile terminating in a smooth, polished cupped blade tip engineered to cradle nerve roots securely without slippage or traction injury. Its ergonomic flat-surfaced handle frame delivers exceptional tactile sensitivity and steady manual control throughout delicate sequestrectomy and canal decompression.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Love Nerve Root Retractor to withstand repetitive hospital sterilization cycles, providing neurosurgical and orthopedic spine centers worldwide with unyielding structural rigidity, high corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)