



# RETRACTOR MAKER

## PRODUCT BROCHURE OF

Alar Retractor

### CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



## Alar Retractor

**SKU** RM-AR-1

**\$29 USD**

### AVAILABLE SIZES & PRICING

| SIZE / OPTION           | SKU / DETAILS                          | PRICE (USD) |
|-------------------------|----------------------------------------|-------------|
| Alar Retractor          | 15 cm overall length, 22 × 13 mm blade | \$29        |
| Protective Sleeve / Box | Individual transit packaging           | \$29        |

### PRODUCT DESCRIPTION

The Alar Retractor is a precision handheld nasal retractor engineered specifically to elevate and support the alar margin and nasal tip during aesthetic and reconstructive rhinoplasty. Designed with a curved distal blade that conforms securely to the nostril border, this instrument provides reliable anatomical exposure without causing ischemic edge trauma or cartilage distortion.

Optimized for delicate facial procedures, this dedicated Alar Retractor configuration features an overall working length of 15 cm with a precisely curved 22 mm wide by 13 mm long blade. Its balanced ergonomic flat handle offers tactile feedback and steady grip control, allowing surgeons and assistant staff to maintain steady upward retraction in deep nasal vestibule corridors without visual obstruction.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Alar Retractor to strict clinical standards, providing cosmetic clinics and ENT hospitals with a fully autoclavable, corrosion-resistant instrument ready for daily operative use.



# THANK YOU

We look forward to discussing your surgical retractor requirements.

## CONTACT RETRACTOR MAKER

EMAIL

**[info@retractor maker.com](mailto:info@retractor maker.com)**

WHATSAPP

**+92 327 8345101**

WEBSITES

**[retractor maker.com](http://retractor maker.com)**  
**[retractor manufacturer.com](http://retractor manufacturer.com)**

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

## FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)