



RETRACTOR MAKER

PRODUCT BROCHURE OF

Hess Nerve Root Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Hess Nerve Root Retractor

SKU RM-HNRR-1

\$69 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
10.0 cm — 3.0 mm - 4.0 mm blade — Standard	RM-HNRR-1-1	\$69
10.0 cm — 5.0 mm - 7.0 mm blade — Standard	RM-HNRR-1-2	\$69
12.0 cm	RM-HNRR-1-3	\$74
17.5 cm	RM-HNRR-1-4	\$79

PRODUCT DESCRIPTION

The Hess Nerve Root Retractor is a specialized neurosurgical and spinal retractor engineered specifically to provide safe, non-traumatic retraction and protection of traversing nerve roots and the thecal sac during lumbar disc procedures. Designed with an angled working end and a cupped tip, this retractor allows spine surgeons to gently mobilize nerve fibers medially, exposing the herniated intervertebral disc without risking traction neuropraxia.

Built to adapt to varying spinal depths, this Hess Nerve Root Retractor series features three distinct physical reaches: 10.0 cm (\$69), 12.0 cm (\$74), and 17.5 cm (\$79). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument offers four selectable blade width options: 3.0 mm, 4.0 mm, 5.0 mm, and 7.0 mm. Its blade structural geometry combines an angled blade architecture with a specialized grip handle assembly that gives operating surgeons exceptional tactile feedback and fatigue-free manual leverage under microscopic illumination.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Hess Nerve Root Retractor to withstand repetitive hospital autoclave cycles, providing neurosurgical suites and orthopedic spine centers

worldwide with unyielding structural rigidity, high corrosion resistance, and complete clinical durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)