



RETRACTOR MAKER

PRODUCT BROCHURE OF

PERCY AMPUTATION RETRACTOR

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



PERCY AMPUTATION RETRACTOR

SKU RM-PAR-1

\$399 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
24.3 cm	RM-PAR-1-1	\$399
35.0 cm	RM-PAR-1-2	\$449
42.0 cm	RM-PAR-1-3	\$499

PRODUCT DESCRIPTION

The PERCY AMPUTATION RETRACTOR is a classic orthopedic and trauma surgical retractor engineered specifically to provide broad, circumferential soft-tissue and muscular retraction during major limb amputations. Designed with an expansive circular shielding paddle that encircles the bone shaft, this instrument allows surgical assistants to pull transected muscle bellies proximal to the bone cut cleanly, protecting soft tissues from rotating or oscillating amputation saws.

Built to offer immense durability and long-life performance in demanding trauma environments, this PERCY AMPUTATION RETRACTOR series is available across three selectable overall lengths: 24.3 cm (\$399) for upper extremity and pediatric procedures, 35.0 cm (\$449) for transfemoral and transtibial amputations, and an extended 42.0 cm (\$499) model for bariatric or heavy muscular limbs. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument features an asymmetrical detachable handle layout that provides maximum pulling leverage while allowing rapid disassembly for easy cleaning.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the PERCY AMPUTATION RETRACTOR to withstand heavy mechanical traction and repeated hospital autoclave cycles, providing trauma centers, military field hospitals, and orthopedic surgical suites worldwide with unyielding structural reliability and complete clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)