



RETRACTOR MAKER

PRODUCT BROCHURE OF

Ferreira Fiber Optic Face Lift Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Ferreira Fiber Optic Face Lift Retractor

SKU RM-FFOFLR-1

\$99 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Ferreira Fiber Optic Face Lift Retractor	4 ¼" (108 mm) length, 22 mm blade, fiber-optic bundle
1	Universal Light Cable Fitting	Connects to standard OR light cables
1	Protective Sleeve / Box	Standard transit packaging

PRODUCT DESCRIPTION

The Ferreira Fiber Optic Face Lift Retractor is a premier illuminated aesthetic surgical instrument engineered specifically to provide superior facial flap elevation combined with direct, shadowless optical illumination during rhytidectomy. By directing cool light into deep surgical tunnels beneath the skin flap, this instrument illuminates subtle anatomical planes, allowing surgeons to dissect with confidence while protecting the facial nerve branches and superficial vessels.

Designed for targeted facial flap deflection and cavity clearance, this Ferreira Fiber Optic Face Lift Retractor model features an overall tool length of 4 ¼ inches (108 mm, \$99) paired with a 22 mm wide blade. Its gently curved blade face distributes traction evenly across cutaneous tissues, preventing focal compression damage. The integrated fiber-optic bundle incorporates a universal fitting compatible with standard hospital light cables (Storz, Olympus, Wolf, ACMI), ensuring effortless connectivity to operating room illumination sources.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) and high-transmission optical glass fibers with an anti-glare satin matte finish, Retractor Maker produces the Ferreira Fiber Optic Face Lift Retractor to strict clinical tolerances, delivering exceptional durability, optical clarity, and complete autoclave durability for cosmetic and facial plastic surgery centers worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)