



RETRACTOR MAKER

PRODUCT BROCHURE OF

Redundant Nonconductive Wall Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Redundant Nonconductive Wall Retractor

SKU RM-RNWR-7IN-1

\$99 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Medium — 7 in, 9.0×2.0 cm blade	RM-RNWR-7IN-1-1	\$99
Large — 7 in, 9.0×2.5 cm blade	RM-RNWR-7IN-1-2	\$109
X-Large — 7 in, 9.0×3.0 cm blade	RM-RNWR-7IN-1-3	\$119

PRODUCT DESCRIPTION

The Redundant Nonconductive Wall Retractor is a specialized gynecological surgical instrument engineered specifically to prevent redundant lateral vaginal wall tissue from collapsing into the operative corridor during loop electrosurgical excision procedures (LEEP) and laser ablation surgery. Coated with a durable, heat-resistant non-conductive insulation layer, this retractor prevents stray electrical currents and thermal arcing from damaging surrounding healthy vaginal mucosa.

Built in a single standard size configuration featuring specialized protective insulation, this Redundant Nonconductive Wall Retractor model features an overall tool length of 7 inches (17.8 cm, \$99). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the working blade measures 9.0 cm long × 2.0 cm wide, providing broad surface support against lateral vaginal folds. The surface profile finish combines highly polished handling interfaces with an electrosurgical-safe blue nylon insulated working layer that withstands repetitive high-frequency cutting and coagulation without dielectric breakdown. The ergonomic flat handle ensures fatigue-free manual leverage for operating assistants.

Forged to strict clinical standards, Retractor Maker manufactures the Redundant Nonconductive Wall Retractor to withstand repetitive hospital autoclave cycles, providing gynecological surgery suites, colposcopy clinics, and women's health centers worldwide with exceptional electrical safety, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)