



RETRACTOR MAKER

PRODUCT BROCHURE OF

Tedan Phantom CS-1600 Cervical Retractor Set

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Tedan Phantom CS-1600 Cervical Retractor Set

SKU RM-TPCS1600-SET-1

\$2,899 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
10 Pairs (20 Total Blades)	Blunt-Edge Retractor Blades	Ranging from 30 mm to 75 mm in lengths (5 mm increments)
6 Pairs (12 Total Blades)	Micro Blunt-Edge Retractor Blades	Ranging from 35 mm to 60 mm in lengths (5 mm increments)
8 Blades	Long-Toothed Retractor Blades	Ranging from 30 mm to 65 mm in lengths (5 mm increments)
1 Pair (2 Total Blades)	Paddle Blades	Standard multi-directional distractor frame compatibility
2 Handles	Side-Loading Retractor Handles	Ergonomic, twist-release design for rapid modular blade changes
1 Gauge	Color-Coded Depth Gauge	Calibrated for rapid soft tissue depth measurement
1 Assembly	Twist Drill with Driver Handle	Manual starter bit assembly for distraction screws
1 Guide	Dual Drill Guide	For securing exact parallel spacing between vertebral distraction screws
1 Guide	Single Drill Guide	For independent or modular distraction screw placement
1 Bit	Twist Drill Bit	Dedicated replacement bit for distraction screw pilot drilling

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
2 Drivers	Distraction Screw Drivers	Secure-grip manual installation drivers for Caspar-style pins
1 Tool	Caliper Measuring Tool	Precision gauge for confirming disc height and implant gap tolerances
1 Tool	Bone Tamp with Guard	Dual-ended surgical seat tool with 6 mm × 10 mm profile
1 Container	Sterilization Box	Heavy-duty, custom-fitted organizational surgical tray

PRODUCT DESCRIPTION

The Tedan Phantom CS-1600 Cervical Retractor Set is an internationally recognized comprehensive anterior cervical surgical retractor and distractor system engineered specifically to provide safe, non-traumatic soft-tissue displacement and intervertebral body distraction during anterior cervical discectomy and fusion (ACDF). Designed with quick-connect side-loading handles and low-profile blade geometries, this system allows spine surgeons to retract the trachea, esophagus, and carotid sheath cleanly, maintaining an unhindered visual corridor down to the disc space without obstructing fluoroscopic imaging or microscopic sightlines.

Built to clinical framework specifications, this Tedan Phantom CS-1600 Cervical Retractor Set complete system (\$2,899) features an extensive instrumentation suite organized in a heavy-duty custom-fitted sterilization container. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the modular blade inventory includes 10 pairs (20 total) of blunt-edge retractor blades from 30 mm to 75 mm, 6 pairs (12 total) of micro blunt-edge blades from 35 mm to 60 mm, 8 long-toothed blades from 30 mm to 65 mm, and 1 pair of multi-directional paddle blades. The set includes 2 ergonomic twist-release side-loading handles, a color-coded depth gauge, a manual twist drill starter bit, dual and single drill guides, replacement drill bits, 2 distraction screw drivers, a precision caliper measuring tool, and a 6 mm × 10 mm guarded bone tamp.

Forged with an anti-glare satin matte finish and color-coded anodized blade sizing markers, Retractor Maker manufactures the Tedan Phantom CS-1600 Cervical Retractor Set to withstand repetitive hospital autoclave cycles, providing neurosurgical suites, orthopedic spine centers, and surgical trauma theaters worldwide with exceptional structural precision, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)