



RETRACTOR MAKER

PRODUCT BROCHURE OF

Freeman Flap Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Freeman Flap Retractor

SKU RM-FFR-8-1

\$19 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
½" (12 mm) Wide — Curved or Straight	RM-FFR-8-1-1	\$19
1" (24 mm) Wide — Curved or Straight	RM-FFR-8-1-2	\$24

PRODUCT DESCRIPTION

The Freeman Flap Retractor is an indispensable handheld plastic surgery instrument engineered specifically to provide steady, atraumatic flap elevation and subcutaneous tissue retraction during aesthetic and reconstructive procedures. Designed with a smooth blade profile terminating in a gentle retaining lip, this instrument allows surgeons to maintain expansive working pockets without causing tissue tearing or edge necrosis.

Built for complete procedural adaptability, this Freeman Flap Retractor model features an overall tool length of 8 inches (20.3 cm). Surgeons can customize the instrument through selectable blade widths: a ½-inch (12 mm) wide blade (\$19) for delicate, localized flap elevation and a 1-inch (24 mm) wide blade (\$24) for broader cutaneous planes. Both widths are available in selectable Curved or Straight tip profiles to match varying tissue vectors. The balanced flat handle provides a comfortable, fatigue-free manual grip.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Freeman Flap Retractor to withstand repetitive hospital autoclave cycles, providing plastic, aesthetic, and general surgery suites worldwide with exceptional corrosion resistance, structural rigidity, and complete clinical reliability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)