



RETRACTOR MAKER

PRODUCT BROCHURE OF

Lobster General Retractor System

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Lobster General Retractor System

SKU RM-LGRS-1

\$2,699 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Lobster Frame	—
1	Lobster Table Clamp with Vertical Post & Frame	—
6	Lobster Retractor Clamp	—
2	Lobster Mayo Swivel Retractor	$\frac{3}{4}$ " $\times$ 2"
1	Lobster Fence Swivel Retractor	4" $\times$ 6"
1	Lobster Splanchnic Swivel Retractor	2" $\times$ 6"
1	Lobster Deaver Swivel Retractor	$2\frac{3}{4}$ " $\times$ $4\frac{1}{2}$ "
1	Lobster Malleable Swivel Retractor (Perforated)	2" $\times$ $6\frac{1}{2}$ "
1	Lobster Malleable Swivel Retractor (Perforated)	3" $\times$ $6\frac{1}{2}$ "
1	Lobster Harrington Sweetheart Retractor	2" $\times$ 5"
1	Lobster Kelly Swivel Retractor, Modified	3" $\times$ $2\frac{1}{2}$ "
1	Lobster Kelly Swivel Retractor, Modified	3" $\times$ 4"
1	Sterilization Box	—

PRODUCT DESCRIPTION

The Lobster General Retractor System is an adaptable table-mounted surgical retractor platform engineered to eliminate assistant fatigue and establish continuous, hands-free exposure across open abdominal procedures. By clamping securely to the operating table side rail, this self-retaining frame creates an unyielding structural perimeter that maintains stable tissue retraction without slipping during lengthy operations.

Configured specifically for general abdominal and visceral interventions, the Lobster General Retractor System includes a heavy-duty table post clamp, a rigid Lobster frame, and six quick-locking slide clamps. The system features a 10-piece general surgery swivel blade assortment—including modified Kelly swivel blades (3" × 2½" and 3" × 4"), deep Splanchnic visceral blades, broad Fence retractors, Deaver, Harrington Sweetheart, Mayo, and perforated malleable strips that flex to patient contours.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with an anti-glare satin finish, Retractor Maker supplies the Lobster General Retractor System non-sterile and fully autoclavable in a custom heavy-duty sterilization box, providing hospitals and surgical centers with reliable durability and consistent OR performance.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)