



RETRACTOR MAKER

PRODUCT BROCHURE OF

Williger Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Williger Retractor

SKU RM-WLR-13-1

\$79 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Williger Retractor (3 Prongs) — 13 cm	RM-WLR-13-1-1	\$79
Williger Retractor (4 Prongs) — 13 cm	RM-WLR-13-1-2	\$89

PRODUCT DESCRIPTION

The Williger Retractor is an essential oral, maxillofacial, and periodontal surgical retractor engineered specifically to provide secure, non-slip soft-tissue and mucoperiosteal flap displacement during bone procedures. Designed with a curved rake head armed with sharp pointed prongs, this instrument allows surgeons to grasp fibrous periosteum and gingival tissues firmly, exposing cortical bone margins without slipping into adjacent vascular structures.

Built for refined tactile control in confined intraoral fields, this Williger Retractor maintains a standardized overall length of 13 cm. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument offers two distinct model configurations: a 3-prong model (\$79) that delivers a highly focused, strong, and controlled grasp of delicate soft tissue margins, and an expanded 4-prong model (\$89) that provides an expanded tight gripping area to ensure maximum exposure of the bone layout. The ergonomic flat handle ensures operating assistants maintain steady, fatigue-free manual leverage throughout demanding bone cuts.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Williger Retractor to withstand repetitive hospital sterilization cycles, providing oral surgery clinics, dental implantology suites, and facial reconstructive theaters worldwide with exceptional corrosion resistance, structural rigidity, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)