



RETRACTOR MAKER

PRODUCT BROCHURE OF

Rollet Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Rollet Retractor

SKU RM-RLR-5IN-1

\$39 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Standard — 13 cm, 4 sharp prongs, 4.5 mm tip	RM-RLR-5IN-1-1	\$39
Large — 13.5 cm, 4 prongs, 7.0 mm tip	RM-RLR-5IN-1-2	\$44

PRODUCT DESCRIPTION

The Rollet Retractor is a specialized handheld surgical rake retractor engineered specifically to provide precise, non-slip skin flap elevation and soft-tissue retraction during delicate ophthalmic, facial plastic, and reconstructive procedures. Designed with finely sharpened hooked prongs, this instrument allows surgeons to grasp cutaneous and subcutaneous layers securely without slipping or crushing fragile marginal tissues.

Classified under operating-room grade hooks and handheld retractors, this Rollet Retractor configuration features an overall tool length of 5 inches (13 cm, \$39). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument is built to withstand the rigors of frequent autoclave sterilization while offering a lifetime warranty against rust. The prong profile configuration incorporates a 4 sharp prongs layout optimized for slippage-free tissue engagement in dense fascial planes and periosteum. Its slender ergonomic handle ensures fine fingertip control.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Rollet Retractor to withstand repetitive hospital sterilization cycles, providing ophthalmic, plastic, and reconstructive surgical suites worldwide with exceptional corrosion resistance, structural rigidity, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)