



RETRACTOR MAKER

PRODUCT BROCHURE OF

Patella Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Patella Retractor

SKU RM-PTR-9IN-1

\$84 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Patella Retractor	9" (22.86 cm) length, 22 mm pad, 16 mm blade, hook tip	\$84
Protective Sleeve / Box	Individual transit packaging	\$84

PRODUCT DESCRIPTION

The Patella Retractor is an indispensable orthopedic surgical instrument engineered specifically to provide safe, non-slip eversion and displacement of the patella and patellar tendon during total knee arthroplasty and joint reconstructions. Designed to hook under the patellar margin, this retractor allows surgical teams to mobilize the extensor mechanism cleanly away from the distal femur, establishing an expansive visual corridor to the anterior cruciate ligament and femoral condyles.

Built to ensure absolute durability and corrosion resistance, this dedicated Patella Retractor configuration features an overall tool length of 9 inches (22.86 cm, \$84). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), its sleek, elongated structural geometry incorporates a 22 mm wide patella pad paired with a 16 mm lower blade width featuring a slightly curved edge and a distinctive hook-like functional tip layout. This anatomical profile seats securely behind cortical bone and soft tissue without crushing patellar articular cartilage. The ergonomic handle delivers fatigue-free manual leverage for operating assistants.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Patella Retractor to withstand heavy orthopedic torque and repeated hospital autoclave cycles, providing sports medicine specialists and knee reconstructive surgeons worldwide with exceptional mechanical strength, corrosion resistance, and complete clinical durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)