



RETRACTOR MAKER

PRODUCT BROCHURE OF

Ambidextrous Self-Retaining Retractors

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Ronald Cameron



Ambidextrous Self-Retaining Retractors

SKU RM-ASRR-AMB-1

\$239 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
4×3 teeth layout	RM-ASRR-AMB-1-1	\$239
4×4 teeth layout	RM-ASRR-AMB-1-2	\$249
7×7 teeth layout	RM-ASRR-AMB-1-3	\$259

PRODUCT DESCRIPTION

The Ambidextrous Self-Retaining Retractors provide an adaptable, dual-sided surgical exposure platform engineered to accommodate both left-handed and right-handed surgical approaches. Built with a symmetrical finger-ring chassis and reversible ratchet catch, these retractors eliminate positioning awkwardness and allow scrub teams to mount the instrument from either surgical table margin.

Designed for robust tissue engagement across shallow and intermediate wound depths, the Ambidextrous Self-Retaining Retractors series is offered across three distinct jaw and prong configurations: a 4×3 teeth layout (\$239), a balanced 4×4 teeth layout (\$249), and a wide 7×7 teeth layout (\$259) designed to distribute holding pressure across broad muscle planes. The multi-step mechanical ratchet locks securely under high tissue tension, preventing accidental release during active dissection.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin finish, Retractor Maker manufactures Ambidextrous Self-Retaining Retractors to exacting clinical standards, ensuring smooth ratchet action, high corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)