



RETRACTOR MAKER

PRODUCT BROCHURE OF

Ramus Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Ramus Retractor

SKU RM-RR-MAX-1

\$69 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Small — 17 cm, 11/15×70 mm blade	RM-RR-MAX-1-1	\$59
Standard — 22 cm, 11×70 mm blade	RM-RR-MAX-1-2	\$69
Large — 22 cm, 12/22×70 mm blade	RM-RR-MAX-1-3	\$79

PRODUCT DESCRIPTION

The Ramus Retractor is an indispensable maxillofacial surgical retractor engineered specifically to provide safe, non-traumatic soft-tissue displacement and osseous exposure of the ascending mandibular ramus during orthognathic interventions. Designed with an anatomical channel curve that wraps over the anterior and lateral ramus borders, this retractor shields the facial nerve branches, masseter muscle, and facial vessels from oscillating saws and rotary burs.

Built for deep intraoral clearance, this Ramus Retractor configuration (\$69) features an elongated shank that keeps the assistant's hands safely outside the surgical corridor, preserving an open visual field for the operating surgeon. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the working end incorporates a contoured channel blade that seats securely against cortical bone without slipping into surrounding soft tissues. The ergonomic handle delivers steady, fatigue-free manual leverage throughout complex mandibular osteotomies.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Ramus Retractor to withstand repetitive hospital sterilization cycles, providing oral, maxillofacial, and facial reconstructive surgical teams worldwide with exceptional structural rigidity, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)