



RETRACTOR MAKER

PRODUCT BROCHURE OF

Valve XS Retractor System

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Valve XS Retractor System

SKU RM-VXRS-1

\$999 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Valve XS Retraction System for Lateral Approach	Double Joint configuration
1	Blade Only with Ball Snap Closure	53 × 36 mm
1	Accessories for Lateral Rib Retractor (Malleable)	60 × 35 mm
1	Valve XS Connecting Clamp Type B	Structural clamp
1	Retraction Rod	9" (230 mm)
1	Retractor for MICS CABG	Dedicated component
1	Applying Instrument for Atrium Lift Retractor Blade	Dedicated applying tool
1	Blade Only	30 × 35 mm
1	Blade Only	40 × 35 mm
1	Blade Only	50 × 35 mm
1	Blade Only	60 × 35 mm
1	Blade Only	70 × 35 mm
1	Sterilization Box	Storage container

PRODUCT DESCRIPTION

The Valve XS Retractor System is an advanced thoracic and cardiac retraction platform engineered specifically for minimally invasive valve procedures performed through a small lateral minithoracotomy. By utilizing a double-jointed articulated spreader frame, this system creates an expansive operative corridor while minimizing intercostal rib spread and preventing rib fractures.

Designed for precision coronary and valve interventions, the Valve XS Retractor System includes a double-joint lateral approach retractor frame, a Type B connecting clamp, a 9-inch (230 mm) retraction rod, a dedicated MICS CABG retractor, and a specialized atrium lift applying instrument. The blade set features a 53 × 36 mm ball-snap blade, a malleable lateral rib accessory (60 × 35 mm), and five progressive blades ranging from 30 × 35 mm to 70 × 35 mm to match patient-specific thoracic depths.

Manufactured entirely from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a fine satin finish, Retractor Maker supplies the Valve XS Retractor System non-sterile and fully autoclavable inside a custom heavy-duty aluminum sterilization container, ensuring rapid tray turnover and long-term durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)