



RETRACTOR MAKER

PRODUCT BROCHURE OF

Brunner Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Brunner Retractor

SKU RM-BRUN-1

\$54 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
89 × 19 mm	RM-BRUN-1-1	\$54
101 × 25 mm	RM-BRUN-1-2	\$56
120 × 25 mm	RM-BRUN-1-3	\$58
130 × 29 mm	RM-BRUN-1-4	\$60
140 × 29 mm	RM-BRUN-1-5	\$62
159 × 29 mm	RM-BRUN-1-6	\$64
178 × 29 mm	RM-BRUN-1-7	\$69

PRODUCT DESCRIPTION

The Brunner Retractor is an essential deep-cavity surgical retractor engineered specifically to provide reliable, hands-free soft-tissue and muscular retraction during major open abdominal, thoracic, and pelvic procedures. Designed with an elongated right-angled blade terminating in a pronounced retaining lip, this instrument allows surgical teams to displace thick tissue layers and large muscle masses with minimal manual strain.

Built for complete procedural adaptability, this Brunner Retractor series features a uniform 25.5 cm overall length paired with seven progressive blade dimensions: 89 × 19 mm (\$54), 101 × 25 mm (\$56), 120 × 25 mm (\$58), 130 × 29 mm (\$60), 140 × 29 mm (\$62), 159 × 29 mm (\$64), and an extra-deep 178 × 29 mm (\$69) configuration. The broad blade face provides uniform pressure distribution to prevent tissue necrosis, while the ergonomic teardrop-ring handle ensures a firm, non-slip grip during prolonged counter-traction.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with an anti-glare satin matte finish, Retractor Maker manufactures the Brunner Retractor to withstand repetitive hospital sterilization cycles, providing general, orthopedic, and gynecological surgical teams with exceptional durability and clinical performance.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)