



RETRACTOR MAKER

PRODUCT BROCHURE OF

Beard Distal Radius Wide Hohmann Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Beard Distal Radius Wide Hohmann Retractor

SKU RM-BDRWHR-1

\$79 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Beard Distal Radius Wide Hohmann Retractor	5.37" (13.63 cm) overall length, wide Hohmann blade	\$79
Protective Sleeve / Box	Individual transit packaging	\$79

PRODUCT DESCRIPTION

The Beard Distal Radius Wide Hohmann Retractor is a precision-engineered orthopedic hand and wrist instrument designed to provide superior bone exposure during distal radius fracture fixation and osteotomy procedures. Designed with an angled shaft that integrates a wide, precision-milled blade assembly, this retractor allows surgical teams to leverage deep soft-tissue layers away from the volar and dorsal radius with minimal manual strain.

Built for optimal wrist surgery ergonomics, this Beard Distal Radius Wide Hohmann Retractor configuration features a compact overall length of 5.37 inches (13.63 cm). Its wide blade profile is contoured to seat stably against the radial shaft, lifting the pronator quadratus and flexor tendons cleanly away from the plating zone. The specialized ergonomic grip handle provides maximum manual control and a secure hold even under fluid exposure, significantly reducing assistant hand fatigue during complex reduction maneuvers.

Manufactured from high-tensile, Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a glare-absorbing satin matte finish, Retractor Maker supplies the Beard Distal Radius Wide Hohmann Retractor non-sterile and fully autoclavable inside protective packaging, delivering high durability and corrosion resistance for orthopedic and hand surgery theaters worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)