



RETRACTOR MAKER

PRODUCT BROCHURE OF

Deep Gelpi Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Deep Gelpi Retractor

SKU RM-DGR-SERIES-1

\$68 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
3 ½" - 4 ½" — Sharp, Blunt, Angled Sharp, 90°, Right Angled — Superficial minor surgery, hand, and pediatric exposure	RM-DGR-SERIES-1-1	\$68
5 ½" - 6 ¾" — Sharp, Blunt, Angled Sharp, 90°, Right Angled — Intermediate soft tissue, plastic, and cervical access	RM-DGR-SERIES-1-2	\$74
7 ½" - 7 ¾" — Sharp, Blunt, Angled Sharp, 90°, Right Angled — Standard lumbar spine laminectomy and joint retraction	RM-DGR-SERIES-1-3	\$84
9 ¼" - 12" — Sharp, Blunt, Angled Sharp, 90°, Right Angled — Deep cavity pelvic, bariatric, and multi-level spinal fusion	RM-DGR-SERIES-1-4	\$94

PRODUCT DESCRIPTION

The Deep Gelpi Retractor is an essential self-retaining surgical instrument engineered specifically to provide reliable, multi-depth tissue retraction across a wide variety of surgical specialties. Built with out-turned prongs that anchor firmly into wound margins, this retractor allows surgical teams to achieve expansive operative access while eliminating manual assistant fatigue.

Designed to accommodate every anatomical depth and clinical scenario, the Deep Gelpi Retractor series features an extensive matrix of customizable options. Surgeons can choose across ten overall lengths: 3 ½" (\$68), 4 ½" (\$72), 5 ½" (\$76), 6 ½" (\$74), 6 ¾" (\$80), 7 ½" (\$84), 7 ¾" (\$88), 9 ¼" (\$94), 11" (\$102), and 12" (\$109). The working tips and arms are selectable in Sharp Points, Blunt Points, Angled Sharp Points, Angled 90

Degree, and Right Angled Arms. Its precision finger-ring ratchet catch locks firmly under tissue counter-tension.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Deep Gelpi Retractor to withstand repetitive hospital autoclave cycles, providing operating rooms worldwide with exceptional mechanical rigidity, smooth ratchet action, and lifelong durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)