



RETRACTOR MAKER

PRODUCT BROCHURE OF

Oldberg retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Oldberg retractor

SKU RM-OR-21-1

\$49 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
13.0 mm Wide Blade	RM-OR-21-1-1	\$49
15.0 mm Wide Blade	RM-OR-21-1-2	\$49
19.0 mm Wide Blade	RM-OR-21-1-3	\$49

PRODUCT DESCRIPTION

The Oldberg retractor is a specialized handheld surgical instrument engineered specifically to provide reliable, non-traumatic soft-tissue displacement and wound margin exposure during cranial, spinal, and deep orthopedic procedures. Designed with an elongated slender shank, this retractor allows surgical assistants to maintain steady manual counter-traction without obstructing the surgeon's direct visual corridor or crowding surgical loupes and operating microscopes.

Built to accommodate varying surgical fields, this Oldberg retractor maintains a uniform longitudinal profile of 21 cm (\$49). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument offers three unique blade footprint widths selectable from the option parameters: a 13.0 mm wide blade featuring a 3-prong rake profile for focused fascial purchase, a 15.0 mm wide blade featuring a 4-prong rake profile for intermediate wound margin displacement, and an expansive 19.0 mm wide blade featuring a smooth, plain blade profile for broad, atraumatic soft-tissue holding. Its ergonomic flat handle delivers fatigue-free manual leverage.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Oldberg retractor to withstand repetitive hospital autoclave cycles, delivering high corrosion resistance, unyielding structural rigidity, and complete autoclave durability for neurosurgical and spine centers worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)