



RETRACTOR MAKER

PRODUCT BROCHURE OF

Hargis Coronoid Retractor With Chain

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Hargis Coronoid Retractor With Chain

SKU RM-HCRWC-1

\$49 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Hargis Coronoid Retractor	19 cm overall length, notched coronoid blade — 1
1	Integrated Retraction Chain	20 cm heavy-duty stainless steel link chain — 1
1	Protective Sleeve / Box	Individual transit packaging — 1

PRODUCT DESCRIPTION

The Hargis Coronoid Retractor With Chain is a specialized maxillofacial surgical retractor engineered specifically to provide reliable, hands-free exposure of the mandibular coronoid process and ascending ramus during orthognathic interventions. Designed with an anatomical notched beak that hooks firmly over the anterior border of the coronoid, this instrument displaces the temporalis muscle tendon safely away from sagittal split cuts and reciprocating bone saws.

Built to endure rigorous operative forces, this Hargis Coronoid Retractor With Chain configuration features an overall tool length of 19 cm (\$49) paired with a heavy-duty 20 cm retraction chain. The chain assembly allows surgical teams to attach standard surgical weights or drape the retractor to establish continuous gravitational counter-traction, eliminating assistant hand fatigue and keeping the visual field unencumbered. The curved blade contours wrap cleanly around the anterior ramus border without lacerating deep oral mucosa.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Hargis Coronoid Retractor With Chain to withstand repetitive hospital sterilization cycles, providing oral, maxillofacial, and craniofacial surgeons worldwide with exceptional corrosion resistance, mechanical stability, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)