



RETRACTOR MAKER

PRODUCT BROCHURE OF

Redler Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Redler Retractor

SKU RM-RDR-5IN-1

\$49 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Redler Retractor	5" (12.7 cm) overall length, tapered functional blade, slender ergonomic handle	\$49
Protective Sleeve / Box	Individual transit packaging	\$49

PRODUCT DESCRIPTION

The Redler Retractor is an adaptable handheld surgical instrument engineered specifically to provide safe, non-traumatic soft-tissue displacement and wound margin retraction during minor general, plastic, and orthopedic surgical procedures. Designed with an elongated slender shank terminating in a tapered blade end, this retractor allows surgical teams to mobilize cutaneous margins and subcutaneous fascial layers cleanly without causing localized crush trauma.

Built in a single standard size configuration engineered for exceptional durability and corrosion resistance, this Redler Retractor configuration features an overall tool length of 5 inches (12.7 cm, \$49). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the working blade profile incorporates a solid, robust metal chassis with a tapered functional end that fits into tight surgical corridors with ease. The handle assembly features a slender, thin, and elongated handle with a sturdy ergonomic grip designed to reduce hand fatigue during extended delicate dissections.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Redler Retractor to withstand repetitive hospital sterilization cycles, delivering high corrosion resistance, unyielding structural rigidity, and complete autoclave durability for outpatient clinics, plastic surgery suites, and emergency surgical centers worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)