



RETRACTOR MAKER

PRODUCT BROCHURE OF

Super-Slide Laminectomy Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Super-Slide Laminectomy Retractor

SKU RM-SSLR-1

\$2,499 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Retractor Frame	Main structural frame of the retractor system
4	Small Blades	2" deep smaller blade options for surgical retraction
4	Medium Blades	3" deep medium-depth blade options
4	Large Blades	4" deep larger-depth blade options
4	Supplementary Blades	4¾" deep extended-depth blade options
1	Sterilization Box	Container for organizing and storing the complete set

PRODUCT DESCRIPTION

The Super-Slide Laminectomy Retractor is a robust, self-retaining surgical retractor engineered to provide powerful, controlled tissue displacement during extensive posterior laminectomy and spinal fusion operations. Built to handle heavy paraspinal musculature without jumping teeth or flexing under counter-traction, this system establishes an expansive, stationary visual corridor over multiple vertebral segments.

Engineered for seamless intraoperative adjustment, this complete Super-Slide Laminectomy Retractor system features a precision-machined sliding ratchet frame that glides effortlessly along its track to apply precise mechanical retraction. The set includes an 18-piece inventory containing four small blades (2" deep), four medium blades (3" deep), four large blades (4" deep), and four supplementary blades (4¾" deep) designed to seat securely against the laminar borders and facet joints.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with an anti-glare satin finish, Retractor Maker supplies the Super-Slide Laminectomy Retractor non-sterile and fully autoclavable inside a

dedicated heavy-duty sterilization box, providing spine surgeons and hospitals with dependable clinical durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)