



RETRACTOR MAKER

PRODUCT BROCHURE OF

McCarthy Mini-Sternotomy Retractor with Universal Lift

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



McCarthy Mini-Sternotomy Retractor with Universal Lift

SKU RM-MMSR-1

\$699 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Chest Spreader Dual Bars	Primary 5-inch modified spreading frame for sternal separation
1	Split Universal Clamp	Specialized mounting clamp for securing frame attachments
1	Universal Clamp	Multi-purpose clamping block for component integration
1	Main Bracket	Primary mounting bracket for supporting auxiliary retractor arms
1	Tube Holder	Secures surgical tubing within the operative field layout
1	Right Angle Transverse Bar	Structural crossbar for positioning auxiliary retractor elements
1	Large Wire Basket Retractor (Malleable)	Broad malleable wire basket for gentle tissue displacement
2	Small Wire Basket Retractor (Malleable)	Compact malleable wire baskets for targeted tissue retraction

PRODUCT DESCRIPTION

The McCarthy Mini-Sternotomy Retractor with Universal Lift is an advanced cardiac retractor system designed specifically for minimally invasive valve and aortic procedures where minimal surgical trauma is required. By combining reliable lateral sternal spreading with independent vertical lifting capabilities, this system elevates the left sternal margin to provide direct visual alignment with the atrium and heart valves.

Engineered for precision thoracic positioning, the McCarthy Mini-Sternotomy Retractor with Universal Lift features a 5-inch modified dual-bar chest spreader frame, an integrated removable Universal Lift assembly, a split universal clamp, a standard universal clamp, a main support bracket, and an auxiliary right-angle transverse bar. The system is equipped with a suction/irrigation tube holder and three malleable wire basket retractors (one large and two small) that adapt to cardiac tissue margins without focal trauma.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a non-reflective satin finish, Retractor Maker supplies the McCarthy Mini-Sternotomy Retractor with Universal Lift non-sterile and fully autoclavable inside a heavy-duty protective storage container, delivering surgical durability and dependable OR performance.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)