



RETRACTOR MAKER

PRODUCT BROCHURE OF

Volkman Retractor-USA Pattern

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Volkman Retractor-USA Pattern

SKU RM-VRUP-223-1

\$64 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
2 Prongs (Blunt / Sharp) — 8 mm - 11 mm	RM-VRUP-223-1-1	\$64
3 Prongs (Blunt / Sharp) — 11 mm - 18 mm	RM-VRUP-223-1-2	\$69
4 Prongs (Blunt / Sharp) — 18 mm - 26 mm	RM-VRUP-223-1-3	\$74
6 Prongs (Blunt / Sharp) — 26 mm - 40 mm	RM-VRUP-223-1-4	\$79

PRODUCT DESCRIPTION

The Volkman Retractor-USA Pattern is a foundational handheld surgical rake retractor engineered specifically to provide reliable, non-slip skin flap elevation and soft-tissue displacement across general, orthopedic, and plastic procedures. Designed with multiple curved rake prongs, this instrument allows surgical assistants to apply uniform counter-traction across cutaneous and subcutaneous planes without causing focal tissue tears or margin ischemia.

Built to withstand high-frequency autoclave sterilization cycles without degrading, this Volkman Retractor-USA Pattern framework maintains a uniform total length of 22.3 cm that gives operating assistants excellent reaching clearance. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument offers available prong quantity footprints with 2 (\$64), 3 (\$69), 4 (\$74), or 6 (\$79) prongs. Surgeons can select between Blunt prong tips engineered for safe, atraumatic soft tissue manipulation, or Sharp prong tips engineered to deliver a firm, slippage-free grip on tougher tissue boundaries. The selectable blade width dimensions span from 8 mm to 26 mm wide, with specific options including 11 mm, 18 mm, 26 mm, and 40 mm.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Volkman Retractor-USA Pattern to uphold strict clinical standards, providing trauma centers, plastic surgery suites, and general operating theaters worldwide with exceptional corrosion resistance, structural rigidity, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



LinkedIn



Facebook



Instagram



X (Twitter)



YouTube



Pinterest