



RETRACTOR MAKER

PRODUCT BROCHURE OF

Mayfield Headrest

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Mayfield Headrest

SKU RM-MHR-1

\$699 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Adult Model	RM-MHR-1-A	\$699
Pediatric Model	RM-MHR-1-P	\$799

PRODUCT DESCRIPTION

The Mayfield Headrest is an essential non-pinned cranial positioning device engineered to deliver rigid head stabilization and superior pressure redistribution during procedures where skeletal skull pins are contraindicated. Designed to connect directly to standard skull clamp base units and operating table adapters, this platform provides reliable immobilization for patients undergoing extensive neurosurgical or spinal interventions.

Engineered to prevent tissue necrosis and localized pressure sores during multi-hour operations, the Mayfield Headrest features a high-density, precision-machined framework paired with fluid-resistant medical-grade gel pads. The platform is offered in two dedicated anatomical options: an Adult horseshoe model featuring independent lateral adjustments and a compact Pediatric configuration specifically contoured to support delicate pediatric cranial structures.

Built to rigorous hospital hygiene standards, Retractor Maker manufactures the Mayfield Headrest with smooth, chemical-resistant surfaces compatible with standard hospital disinfectants, shipping each unit fully tested inside an export-grade protective box with custom-contoured high-density foam packaging.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)