



RETRACTOR MAKER

PRODUCT BROCHURE OF

SHADOWLINE Trimline Cervical Spinal Retractor Set

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



SHADOWLINE Trimline Cervical Spinal Retractor Set

SKU RM-STCSRS-1

\$2,999 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Color-Coded Depth Gauge	Instruments & Frames — standard
2	Locking Retractor Handle (Side-loading)	Instruments & Frames — standard
1	Step-Down Transverse Retractor Frame	Instruments & Frames — double-hinged, side-loading
1	Transverse Retractor Frame	Instruments & Frames — double-hinged, side-loading (120 mm opening)
1	Instrument Case	Storage & Protection — organizing tray
1	Sterilization Box	Storage & Protection — perforated container
2	Silver Tissue Blade	Blunt tissue blades (pairs) — 30 mm × 23 mm
2	Blue Tissue Blade	Blunt tissue blades (pairs) — 35 mm × 23 mm
2	Red Tissue Blade	Blunt tissue blades (pairs) — 40 mm × 23 mm
2	Green Tissue Blade	Blunt tissue blades (pairs) — 45 mm × 23 mm
2	Violet Tissue Blade	Blunt tissue blades (pairs) — 50 mm × 23 mm
2	Black Tissue Blade	Blunt tissue blades (pairs) — 55 mm × 23 mm
2	Gold Tissue Blade	Blunt tissue blades (pairs) — 60 mm × 23 mm
2	Navy Tissue Blade	Blunt tissue blades (pairs) — 65 mm × 23 mm

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
2	Light Green Tissue Blade	Blunt tissue blades (pairs) — 70 mm × 23 mm
2	Orange Tissue Blade	Blunt tissue blades (pairs) — 75 mm × 75 mm
2	Blue Narrow Muscle Blade	Narrow muscle blades (short teeth, pairs) — 35 mm × 14 mm
2	Red Narrow Muscle Blade	Narrow muscle blades (short teeth, pairs) — 40 mm × 14 mm
2	Green Narrow Muscle Blade	Narrow muscle blades (short teeth, pairs) — 45 mm × 14 mm
2	Violet Narrow Muscle Blade	Narrow muscle blades (short teeth, pairs) — 50 mm × 14 mm
2	Black Narrow Muscle Blade	Narrow muscle blades (short teeth, pairs) — 55 mm × 14 mm
2	Gold Narrow Muscle Blade	Narrow muscle blades (short teeth, pairs) — 60 mm × 14 mm
1	Silver Lateral Muscle Blade	Lateral muscle blades (short teeth, singles) — 30 mm × 23 mm
1	Blue Lateral Muscle Blade	Lateral muscle blades (short teeth, singles) — 35 mm × 23 mm
1	Red Lateral Muscle Blade	Lateral muscle blades (short teeth, singles) — 40 mm × 23 mm
1	Green Lateral Muscle Blade	Lateral muscle blades (short teeth, singles) — 45 mm × 23 mm
1	Violet Lateral Muscle Blade	Lateral muscle blades (short teeth, singles) — 50 mm × 23 mm
1	Black Lateral Muscle Blade	Lateral muscle blades (short teeth, singles) — 55 mm × 23 mm
1	Gold Lateral Muscle Blade	Lateral muscle blades (short teeth, singles) — 60 mm × 23 mm
1	Navy Lateral Muscle Blade	Lateral muscle blades (short teeth, singles) — 65 mm × 23 mm
1	Silver Median Muscle Blade	Median muscle blades (long teeth, singles) — 30 mm × 23 mm
1	Blue Median Muscle Blade	Median muscle blades (long teeth, singles) — 35 mm × 23 mm

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Red Median Muscle Blade	Median muscle blades (long teeth, singles) — 40 mm × 23 mm
1	Green Median Muscle Blade	Median muscle blades (long teeth, singles) — 45 mm × 23 mm
1	Violet Median Muscle Blade	Median muscle blades (long teeth, singles) — 50 mm × 23 mm
1	Black Median Muscle Blade	Median muscle blades (long teeth, singles) — 55 mm × 23 mm
1	Gold Median Muscle Blade	Median muscle blades (long teeth, singles) — 60 mm × 23 mm
1	Navy Median Muscle Blade	Median muscle blades (long teeth, singles) — 65 mm × 23 mm
2	Blue Narrow Tissue Blade	Narrow tissue blades (blunt, pairs) — 35 mm × 14 mm
2	Red Narrow Tissue Blade	Narrow tissue blades (blunt, pairs) — 40 mm × 14 mm
2	Green Narrow Tissue Blade	Narrow tissue blades (blunt, pairs) — 45 mm × 14 mm
2	Violet Narrow Tissue Blade	Narrow tissue blades (blunt, pairs) — 50 mm × 14 mm
2	Black Narrow Tissue Blade	Narrow tissue blades (blunt, pairs) — 55 mm × 14 mm
2	Gold Narrow Tissue Blade	Narrow tissue blades (blunt, pairs) — 60 mm × 14 mm

PRODUCT DESCRIPTION

The SHADOWLINE Trimline Cervical Spinal Retractor Set is an advanced 67-piece surgical instrument platform designed to meet the rigorous exposure and soft-tissue retraction demands of anterior cervical spine procedures. Engineered with a low-profile architecture and color-coded blade sizing, this system provides unyielding surgical field visualization while minimizing focal pressure on delicate visceral structures like the esophagus and trachea.

Configured for complete procedural adaptability, the SHADOWLINE Trimline Cervical Spinal Retractor Set features two double-hinged transverse retractor frames (including a step-down frame and a 120 mm wide-opening frame), two side-loading locking handles, and an intuitive color-coded depth gauge. The set includes 61 color-coded blades across blunt tissue, narrow tissue, short-tooth lateral muscle, and long-tooth median muscle geometries, allowing surgeons to match patient anatomy with millimeter accuracy.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) and high-strength metal alloys with an anti-glare satin finish, Retractor Maker supplies the complete SHADOWLINE Trimline Cervical Spinal Retractor Set non-sterile and fully autoclavable inside a heavy-duty aluminum sterilization container with dedicated organizing trays, ensuring seamless hospital tray turnover.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)