



RETRACTOR MAKER

PRODUCT BROCHURE OF

Kirschner Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Kirschner Retractor

SKU RM-KR-SET-1

\$449 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Self-Retaining Guide Rail Base Frame	Oval or Square chassis (240×240 mm or 340×240 mm)
4 to 6 Blades	Modular Blunt Tissue Blades	Choice of 2-Blade Pair Set or 4-Blade Full Set
4	Individual Retention Pin Locks	Slide-and-lock anchoring mechanisms to secure blades
2	Ergonomic Retracting Forceps Beams	Ultra-lightweight structured handles with smooth outer edges
1	Protective Storage / Box	Standard transit packaging

PRODUCT DESCRIPTION

The Kirschner Retractor is an internationally acclaimed self-retaining surgical retractor system engineered specifically to provide broad, multi-directional abdominal wall and visceral displacement during major general, gynecological, and urological laparotomies. Designed with an expansive perimeter base frame, this instrument allows surgical teams to position multiple modular blades at custom angles, establishing an unencumbered, static visual corridor into deep peritoneal and retroperitoneal spaces without requiring manual assistant holding.

Built to adapt to diverse patient anatomies and incisional requirements, this Kirschner Retractor series offers selectable frame blueprint profiles in a 240 mm × 240 mm chassis setup or an extended 340 mm × 240 mm chassis setup. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), surgeons can select between a 2-Blade Pair Set (\$449) consisting of two 40 mm × 65 mm blades and two 50 mm × 65 mm blades, or upgrade to the 4-Blade Full Set (\$499) consisting of three 50 mm × 65 mm blades, one 80 mm × 90 mm blade, one 98 mm × 50 mm blade, and one 105 mm × 35 mm blade. Individual retention pin locks anchor

each blade along the perimeter frame to eliminate intraoperative slippage.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Kirschner Retractor to withstand repetitive hospital autoclave cycles, providing laparotomy suites, oncologic surgery centers, and pelvic reconstructive theaters worldwide with exceptional mechanical stability, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



LinkedIn



Facebook



Instagram



X (Twitter)



YouTube



Pinterest