



RETRACTOR MAKER

PRODUCT BROCHURE OF

Deaver Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Deaver Retractor

SKU RM-DR-STD-1

\$59 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
18 cm L × 10 mm / 19 mm W — Pediatric, shallow abdominal, and thyroid exposure	RM-DR-STD-1-1	\$59
20 cm L × 13 mm / 16 mm / 26 mm W — Standard appendectomy and minor laparotomy access	RM-DR-STD-1-2	\$63
21 cm L × 23 mm W — Intermediate abdominal wall and pelvic retraction	RM-DR-STD-1-3	\$69
23 cm L × 26 mm W — Standard adult exploratory laparotomy exposure	RM-DR-STD-1-4	\$71
26 cm L × 26 mm W — Deep pelvic and gynecological access	RM-DR-STD-1-5	\$73
30 cm L × 26 mm / 38 mm / 50 mm W — Broad visceral and liver displacement	RM-DR-STD-1-6	\$75
30 cm L × 76 mm / 101 mm W — Extra-wide abdominal wall and bariatric visceral hold	RM-DR-STD-1-7	\$81
33 cm L × 26 mm W — Deep retroperitoneal and pelvic cancer surgery	RM-DR-STD-1-8	\$85
36 cm L × 26 mm W — Maximum-reach bariatric and aortic exposure	RM-DR-STD-1-9	\$87

PRODUCT DESCRIPTION

The Deaver Retractor is an essential surgical retractor recognized worldwide for providing broad, atraumatic soft-tissue and visceral displacement during major abdominal, thoracic, and pelvic operations. Designed with a long, flat blade that arches into a deep curve, this instrument allows surgical assistants to retract the abdominal wall, liver, spleen, and intestines smoothly away from the operative corridor.

Built for complete procedural adaptability across variable anatomical depths, the Deaver Retractor series is manufactured in 15 distinct size combinations: 18 cm × 10 mm (\$59), 18 cm × 19 mm (\$61), 20 cm × 13 mm (\$63), 20 cm × 16 mm (\$65), 20 cm × 26 mm (\$67), 21 cm × 23 mm (\$69), 23 cm × 26 mm (\$71), 26 cm × 26 mm (\$73), 30 cm × 26 mm (\$75), 30 cm × 38 mm (\$77), 30 cm × 50 mm (\$79), 30 cm × 76 mm (\$81), 30 cm × 101 mm (\$83), 33 cm × 26 mm (\$85), and 36 cm × 26 mm (\$87). The smooth, polished blade edges prevent mucosal and serosal tearing.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Deaver Retractor to withstand heavy mechanical loads and repeated hospital autoclaving cycles, providing operating theaters worldwide with exceptional durability and clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)