



RETRACTOR MAKER

PRODUCT BROCHURE OF

Kolbel Shoulder Retractor Frame

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Kolbel Shoulder Retractor Frame

SKU RM-KSRF-23-1

\$159 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Small — 36 mm x 36 mm blade	RM-KSRF-23-1-1	\$159
Medium — 53 mm x 36 mm blade	RM-KSRF-23-1-2	\$169
Large — 68 mm x 36 mm blade	RM-KSRF-23-1-3	\$179

PRODUCT DESCRIPTION

The Kolbel Shoulder Retractor Frame is an internationally recognized self-retaining orthopedic instrument engineered specifically to provide reliable, multi-directional soft-tissue exposure during total and reverse shoulder arthroplasty. Designed to host interchangeable modular blades, this framework allows surgeons to retract the deltoid, pectoralis major, and conjoined tendon simultaneously, establishing an expansive, hands-free surgical corridor over the glenohumeral joint.

Built as a comprehensive kit to accommodate varying structural depths during retraction, this Kolbel Shoulder Retractor Frame features an overall instrument length of 23 cm (9 inches, \$159). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the chassis is equipped with a system of free-moving pivoting arms that articulate smoothly to conform to patient-specific shoulder anatomy. The integrated mechanical ratchet gears maintain constant counter-traction under heavy muscular tension, eliminating manual assistant fatigue during glenoid preparation.

Forged with an anti-glare satin matte finish, Retractor Maker produces the Kolbel Shoulder Retractor Frame to withstand heavy mechanical loads and repeated hospital autoclave cycles, delivering high corrosion resistance, smooth gear tracking, and complete clinical durability for shoulder reconstruction and sports medicine suites worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)