



RETRACTOR MAKER

PRODUCT BROCHURE OF

Concave Knee Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Concave Knee Retractor

SKU RM-CKR-1

\$105 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Standard Blade	RM-CKR-1-1	\$105
Narrow Long Blade	RM-CKR-1-2	\$105

PRODUCT DESCRIPTION

The Concave Knee Retractor is a specialized orthopedic surgical retractor engineered specifically to provide exceptional bone leverage and soft-tissue protection during total knee replacement and complex knee reconstructions. Designed with an anatomically concave blade profile, this instrument conforms smoothly around the femoral condyles and tibial plateau, allowing surgical teams to retract collateral ligaments and periarticular soft tissues without slipping.

Built for heavy joint leverage, this Concave Knee Retractor model features a uniform 9.5-inch (24.13 cm) overall structural length (\$105) that ensures optimal working depth while keeping assistant hands well outside the saw path. Surgeons can select between two dedicated blade geometries: a Standard Blade for broad soft-tissue displacement or a Narrow Long Blade optimized for deep, confined joint spaces and minimally invasive knee approaches. The ergonomic handle provides superior manual leverage, reducing user fatigue during long bone cuts.

Forged from high-tensile, Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Concave Knee Retractor to endure heavy orthopaedic torque and repeated hospital autoclave cycles, providing knee reconstruction teams with unyielding clinical durability and performance.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)