



RETRACTOR MAKER

PRODUCT BROCHURE OF

Cobra Style Femur Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Cobra Style Femur Retractor

SKU RM-CSFR-1

\$84 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Cobra Style Femur Retractor	11.75" (29.84 cm) length, cobra-curved blade, large hooking tip	\$84
Protective Sleeve / Box	Individual transit packaging	\$84

PRODUCT DESCRIPTION

The Cobra Style Femur Retractor is a specialized orthopedic surgical retractor engineered specifically to provide deep, uncompromised exposure of the femoral shaft and proximal femur during major hip reconstructions and fracture fixations. Built with an anatomical cobra-like curve, this instrument allows surgical assistants to leverage thick quadriceps and hamstring musculature smoothly away from the operative corridor.

Designed for robust deep-tissue displacement, this Cobra Style Femur Retractor configuration features an overall length of 11.75 inches (29.84 cm, \$84). Its working end incorporates a large, smooth hooking tip specifically engineered to seat securely around the posterior and medial femoral cortex without slipping or traumatizing surrounding neurovascular bundles. The broad blade body distributes leverage force evenly, preventing focal muscle crushing during heavy retraction.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Cobra Style Femur Retractor to withstand immense orthopaedic torque and repeated hospital autoclaving cycles, delivering exceptional strength, corrosion resistance, and clinical longevity for trauma and joint replacement theaters.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)