



RETRACTOR MAKER

PRODUCT BROCHURE OF

Langenbeck Green Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Langenbeck Green Retractor

SKU RM-LGR-1

\$44 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
17 cm	RM-LGR-1-1	\$44
and 22 cm	RM-LGR-1-2	\$49

PRODUCT DESCRIPTION

The Langenbeck Green Retractor is an indispensable handheld surgical retractor engineered specifically to provide safe, atraumatic soft-tissue and muscle displacement during thyroid, neck, and general surgical procedures. Combining an elongated right-angled blade with a downward-curved retaining lip, this instrument allows surgical assistants to maintain constant counter-traction on delicate tissue planes without causing local bruising.

Built with an adaptive size matrix, the Langenbeck Green Retractor is offered in two overall lengths: 17 cm (\$44) for superficial neck dissections and 22 cm (\$49) for deeper operative fields. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument provides two matching blade footprints: a narrow 25 mm long × 6 mm wide blade for tight thyroid corridors and an intermediate 41 mm long × 11 mm wide blade for broader muscular elevation. Its flat ergonomic handle ensures steady, fatigue-free manual grip control.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Langenbeck Green Retractor to withstand repetitive hospital autoclave cycles, providing operating rooms, ENT suites, and outpatient surgery centers worldwide with exceptional durability, corrosion resistance, and complete clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)