



# RETRACTOR MAKER

## PRODUCT BROCHURE OF

Anderson-Adson Self-Retaining Cerebellar Retractor

### CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



## Anderson-Adson Self-Retaining Cerebellar Retractor

**SKU** RM-AASCR-1

**\$199 USD**

### AVAILABLE SIZES & PRICING

| SIZE / OPTION                                            | SKU / DETAILS | PRICE (USD) |
|----------------------------------------------------------|---------------|-------------|
| Small — 6.5 in, 4x4 sharp prongs                         | RM-AASCR-1-1  | \$189       |
| Standard — 7.5 in (19 cm), 4x4 sharp prongs, 13 mm depth | RM-AASCR-1-2  | \$199       |
| Large — 9.5 in, 4x4 sharp prongs                         | RM-AASCR-1-3  | \$214       |

### PRODUCT DESCRIPTION

The Anderson-Adson Self-Retaining Cerebellar Retractor is a specialized neurosurgical instrument engineered to provide powerful, hands-free scalp and muscle retraction during posterior cranial and high cervical operations. By exerting continuous, locked counter-traction on dense paraspinal and suboccipital musculature, this instrument frees surgical staff and prevents operative field obstruction.

Built for deep posterior fossa exposure, this Anderson-Adson Self-Retaining Cerebellar Retractor configuration features an overall length of 19 cm and an angled frame geometry that keeps the instrument handles positioned flat against patient contours. The working ends feature 4x4 sharp prongs with an integrated 13.0 mm blade depth profile that embeds securely into thick fascia without slipping. The mechanical finger-ring ratchet mechanism allows fine-tuned spread adjustments that remain locked under heavy tissue tension.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a glare-resistant satin matte finish, Retractor Maker manufactures the Anderson-Adson Self-Retaining Cerebellar Retractor to strict clinical tolerances, ensuring smooth ratchet engagement, high corrosion resistance, and complete autoclave durability.



# THANK YOU

We look forward to discussing your surgical retractor requirements.

## CONTACT RETRACTOR MAKER

EMAIL

**[info@retractormaker.com](mailto:info@retractormaker.com)**

WHATSAPP

**+92 327 8345101**

WEBSITES

**[retractormaker.com](http://retractormaker.com)**  
**[retractormanufacturer.com](http://retractormanufacturer.com)**

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

## FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)