



RETRACTOR MAKER

PRODUCT BROCHURE OF

Forearm Swivel Blade Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Forearm Swivel Blade Retractor

SKU RM-FSBR-1

\$74 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Forearm Swivel Blade Retractor	5.12" (13.00 cm) blade, finger-ring handle, ratchet assembly — 1	\$74
Protective Sleeve / Box	Individual transit packaging — 1	\$74

PRODUCT DESCRIPTION

The Forearm Swivel Blade Retractor is a specialized orthopedic extremity instrument engineered specifically to provide targeted joint, soft-tissue, and bone distraction during complex radius and ulna reconstructions. Designed with pivoting swivel blades, this retractor adapts dynamically to shifting anatomical planes, allowing surgeons to maintain expansive, hands-free operative access without exerting uneven shear forces on surrounding muscle tissue.

Built for targeted upper extremity joint and bone distraction, this Forearm Swivel Blade Retractor configuration features a blade length of 5.12 inches (13.00 cm, \$74). The chassis incorporates an ergonomic finger ring design paired with an integrated multi-tooth ratchet assembly that holds targeted tissue tension automatically, eliminating manual assistant fatigue. The smooth, contoured swivel blades conform to variable osseous contours, preventing accidental slippage during fracture reduction and plate fixation.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Forearm Swivel Blade Retractor to withstand repetitive hospital sterilization cycles, providing orthopedic trauma and hand surgery centers worldwide with exceptional mechanical stability, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)