



RETRACTOR MAKER

PRODUCT BROCHURE OF

Sauerbruch Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Sauerbruch Retractor

SKU RM-SBR-9IN-1

\$69 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
$\frac{5}{8}$ " $\times$ 1 $\frac{3}{4}$ " (approx. 16 $\times$ 44 mm) — 9" (22.86 cm)	RM-SBR-9IN-1-1	\$69
2 $\frac{1}{8}$ " $\times$ $\frac{7}{8}$ " (approx. 54 $\times$ 22 mm) — 9" (22.86 cm)	RM-SBR-9IN-1-2	\$69
$\frac{7}{8}$ " $\times$ 2 $\frac{1}{2}$ " (approx. 22 $\times$ 64 mm) — 9" (22.86 cm)	RM-SBR-9IN-1-3	\$69
$\frac{3}{4}$ " $\times$ 2 $\frac{7}{8}$ " (approx. 19 $\times$ 73 mm) — 9" (22.86 cm)	RM-SBR-9IN-1-4	\$69

PRODUCT DESCRIPTION

The Sauerbruch Retractor is an internationally recognized deep-cavity surgical retractor engineered specifically to provide reliable, atraumatic soft-tissue displacement and muscular retraction during thoracic, abdominal, and orthopedic procedures. Designed with an elongated right-angled blade, this instrument allows surgical assistants to maintain constant counter-traction on dense muscle bellies, fascial layers, and visceral walls without slipping into the operative field.

Built to withstand frequent autoclave sterilization cycles while providing a lifetime warranty against rust and defects, this Sauerbruch Retractor features a fixed longitudinal frame of 9 inches (22.86 cm, \$69). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument offers four distinct blade size dimensions selectable from the option parameter menu: $\frac{5}{8}$ " $\times$ 1 $\frac{3}{4}$ ", 2 $\frac{1}{8}$ " $\times$ $\frac{7}{8}$ ", $\frac{7}{8}$ " $\times$ 2 $\frac{1}{2}$ ", and $\frac{3}{4}$ " $\times$ 2 $\frac{7}{8}$ " wide right-angle blade layouts. The smooth, rounded blade margins prevent localized ischemic pressure trauma, while the ergonomic flat handle provides fatigue-free manual leverage.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Sauerbruch Retractor to withstand heavy mechanical loads and repeated hospital sterilization cycles, delivering high corrosion resistance, unyielding structural rigidity, and complete autoclave durability for general, thoracic, and trauma operating theaters worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)