



RETRACTOR MAKER

PRODUCT BROCHURE OF

Cervical Retractor Frame

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Cervical Retractor Frame

SKU RM-CRF-1

\$109 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
16.5 cm Rigid Arms	RM-CRF-1-1	\$109
25.5 cm Hinged Arms	RM-CRF-1-2	\$119

PRODUCT DESCRIPTION

The Cervical Retractor Frame is an essential self-retaining spinal retractor chassis engineered specifically to provide stable, longitudinal soft-tissue exposure during anterior cervical discectomy and fusion (ACDF). Designed to hold interchangeable blunt and toothed blades securely, this frame spreads longus colli muscles and retracts visceral neck structures away from the disc space without manual holding.

Built to accommodate variable clinical approaches and cervical anatomies, this Cervical Retractor Frame is manufactured in two distinct structural configurations. The 16.5 cm Length Frame (\$109) is formatted explicitly with rigid arms to provide a highly consistent, non-yielding retraction force that minimizes involuntary micro-movement during single-level interventions. The extended 25.5 cm Length Frame (\$119) is engineered with hinged arms that articulate dynamically to adapt to varying, asymmetrical patient contours during multi-level procedures.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Cervical Retractor Frame to withstand repetitive hospital sterilization cycles, providing spine surgeons and neurosurgical teams with unyielding rigidity, smooth ratchet tracking, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)