



RETRACTOR MAKER

PRODUCT BROCHURE OF

Sachs Vein Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Sachs Vein Retractor

SKU RM-SVR-22-1

\$54 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Plain Blade Model — 22 cm length, smooth curved blade	RM-SVR-22-1-1	\$54
Toothed Blade Model — 22 cm length, toothed blade edge	RM-SVR-22-1-2	\$54
Sharp Hook Model — 22 cm length, sharp curved hook	RM-SVR-22-1-3	\$54

PRODUCT DESCRIPTION

The Sachs Vein Retractor is a specialized vascular and neurosurgical retractor engineered specifically to provide safe, non-traumatic retraction and mobilization of delicate veins, nerve trunks, and vascular bundles during intricate surgical dissections. Designed with an elongated slender shank, this instrument allows surgeons to gently hold back fragile venous walls without causing focal endothelial bruising or thrombosis.

Built with a uniform body framework of 22 cm (\$54), this Sachs Vein Retractor configuration offers three distinct working ends selectable via the option menu parameters. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) to provide a lifetime warranty against rust and manufacturing defects, surgeons can choose between a Plain Blade model layout for atraumatic vessel cradling, a Toothed Blade model layout engineered to secure an optimal tissue grip and minimize slippage on fibrous adventitia, or a Sharp Hook model layout for precise point-specific retraction. The handle assembly features an integrated ergonomic finger hold layout that ensures steady, fatigue-free manual leverage.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Sachs Vein Retractor to withstand repetitive hospital autoclave cycles, providing vascular, cardiovascular, and neurosurgical operating theaters worldwide with exceptional durability, corrosion resistance, and complete clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)