



RETRACTOR MAKER

PRODUCT BROCHURE OF

Hasson S Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Hasson S Retractor

SKU RM-HSR-LAP-1

\$24 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Rigid Design — 10.0 mm or 13.0 mm	RM-HSR-LAP-1-1	\$24
Malleable Design — 10.0 mm or 13.0 mm	RM-HSR-LAP-1-2	\$44

PRODUCT DESCRIPTION

The Hasson S Retractor is a specialized handheld surgical instrument engineered specifically to provide safe, atraumatic fascial exposure and abdominal wall elevation during the open cut-down (Hasson) laparoscopic technique. Designed with an S-shaped dual-ended architecture, this retractor allows laparoscopic surgeons to elevate umbilical fascial margins cleanly, providing direct visualization for safe peritoneal entry without risking bowel injury.

Built for clinical versatility, this Hasson S Retractor model maintains a fixed overall instrument length of 18.5 cm. The instrument offers customizable material layout profiles: a Malleable Design (\$44) that allows customized intraoperative shaping to meet unique patient contours, and a Rigid Design (\$24) providing maximum structural resistance against dense tissue layers. Both material profiles are available in 10.0 mm and 13.0 mm blade widths with rounded atraumatic edges that protect delicate peritoneum.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Hasson S Retractor to withstand repetitive hospital autoclave cycles, providing laparoscopic and general surgery suites worldwide with exceptional corrosion resistance, structural rigidity, and complete clinical durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)