



RETRACTOR MAKER

PRODUCT BROCHURE OF

Cronin Cheek Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Cronin Cheek Retractor

SKU RM-CCR-ORAL-1

\$24 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Cronin Cheek Retractor	16 cm overall length, 19 mm wide curved blade	\$24
Protective Sleeve / Box	Individual transit packaging	\$24

PRODUCT DESCRIPTION

The Cronin Cheek Retractor is a specialized oral and maxillofacial surgical retractor engineered specifically to provide safe, non-traumatic displacement of the cheeks, lips, and buccal mucosa during intraoral and reconstructive facial surgery. Designed with a broad, curved working end, this instrument allows surgeons to maintain expansive mechanical exposure while shielding soft tissues from rotating rotary burrs and oscillating bone saws.

Built to optimize mechanical exposure and surgical sightlines, this Cronin Cheek Retractor configuration features an overall length of 16 cm (\$24). Its working blade incorporates a 19 mm wide widened blade design configured to smoothly manage and leverage substantial tissue fields without causing localized mucosal bruising or ischemia. The contoured shank keeps the assistant's hands safely outside the surgical field, while the ergonomic handle provides fatigue-free manual leverage throughout extended oral procedures.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Cronin Cheek Retractor to withstand repetitive hospital sterilization cycles, providing oral, maxillofacial, and dental surgical centers worldwide with superior corrosion resistance, unyielding structural rigidity, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)